


map (IProfilerMap)

b2bOpticJobData → items → item → pair → patient → iProfilerData → result → wavefrontCornea → map



point	
occurs	1..n
x	
type	float
occurs	1
unity	µm
description	
y	
type	float
occurs	1
unity	µm
description	
z	
type	float
occurs	1
unity	µm
description	

```
<xs:complexType name="IProfilerMap">
  <xs:sequence maxOccurs="unbounded">
    <xs:element name="point">
      <xs:complexType>
        <xs:sequence>
          <xs:element name="x" type="xs:float" />
          <xs:element name="y" type="xs:float" />
          <xs:element name="z" type="xs:float" />
        </xs:sequence>
      </xs:complexType>
    </xs:element>
  </xs:sequence>
</xs:complexType>
```

From:
<https://wiki2.b2boptic.com/> - **wiki2.b2bOptic.com**

Permanent link:
<https://wiki2.b2boptic.com/en:jobdata:version010603:complextypes:iprofilermap:wavefrontcornea>

Last update: **2017/06/19 18:31**

