


# autorefractorEye (RXDataTypeSimple)

b2boptic → items → item → pair → patient → iProfilerData → result → autorefractorEye



| sphere      |          |
|-------------|----------|
| type        | float    |
| unity       | dpt      |
| occurs      | 1        |
| description |          |
| cylinder    |          |
| type        | Cylinder |
| occurs      | 0..1     |
| description |          |

```
<xs:complexType name="RXDataTypeSimple">
  <xs:sequence>
    <xs:element name="sphere" type="xs:float" />
    <xs:element minOccurs="0" name="cylinder" type="Cylinder" />
  </xs:sequence>
</xs:complexType>
```

From:

<https://wiki2.b2boptic.com/> - **wiki2.b2bOptic.com**

Permanent link:

<https://wiki2.b2boptic.com/en:lensorder:version010602:complextypes:rxdatatypesimple>

Last update:

**2016/05/04 13:32**