


PachyEye (IProfilerResult)

b2boptic → items → item → pair → patient → iProfilerData → result → pachyeye



distCorneaLens	
type	float
unity	dpt
occurs	1
description	distance Cornea Lens

ISNEWLO010604 centralCorneaThickness	
type	float
unity	mm
occurs	0..1
description	Cornea Thickness in middle

```
<xs:complexType name="PachyEye">
  <xs:sequence>
    <xs:element name="distCorneaLens" type="xs:float"/>
    <xs:element name="centralCorneaThickness" type="xs:float"
minOccurs="0"/>
  </xs:sequence>
</xs:complexType>
```

From:
<https://wiki2.b2boptic.com/> - **wiki2.b2bOptic.com**

Permanent link:
<https://wiki2.b2boptic.com/en:lensorder:version010605:complextypes:pachyeye>

Last update: **2024/07/09 04:10**